

In Vitro Diagnostic Reagent

Coagulation Analyzer Cleaning Solution–SY I

This product is used to clean the reaction system of the automated coagulation analyzer during the testing process, to facilitate the in vitro testing of the substance to be measured



Product Photo

Product Model	Coagulation Analyzer Cleaning Solution–SYI
Packaging Specifications	50 mL, 200 mL, 500 mL, 1 L, 20 L
Product Code	p1061

Intended Use

This product is used to clean the reaction system of the automated coagulation analyzer during the testing process, to facilitate the in vitro testing of the substance to be measured

Principle

The sodium hypochlorite in the cleaning solution can effectively decompose the proteins in the reaction system, making the proteins in the blood sample easier to clean...

This document is compiled from product registration information for reference only. Refer to the product instructions and the registration certificate for details.