

In Vitro Diagnostic Reagent

Coagulation Analyzer Cleaning Solution-ZCWY II

This product is a cleaning solution for the reaction system of the automated coagulation analyzer during the testing process, to facilitate the in vitro testing of the substance to be measured



Product Photo

Product Model	Coagulation Analyzer Cleaning Solution- ZCWY II
Packaging Specifications	60 mL, 100 mL, 500 mL, 20 L
Product Code	p1097

Intended Use

This product is a cleaning solution for the reaction system of the automated coagulation analyzer during the testing process, to facilitate the in vitro testing of the substance to be measured

Principle

This product is a cleaning solution for the reaction system. Its main working principle is that the sodium hypochlorite in the cleaning solution can effectively dissolve the proteins, thereby making the proteins in blood specimens easy to clean...

This document is compiled from product registration information for reference only. Refer to the product instructions and the registration certificate for details.