

## In Vitro Diagnostic Reagent

## Coagulation Analyzer Cleaning Solution-ZCWY I

This product is a cleaning solution for the reaction system of the automated coagulation analyzer during the testing process, to facilitate the in vitro testing of the substance to be measured



Product Photo

|               |                                                      |
|---------------|------------------------------------------------------|
| Product Model | Coagulation Analyzer<br>Cleaning Solution-<br>ZCWY I |
|---------------|------------------------------------------------------|

|                             |                   |
|-----------------------------|-------------------|
| Packaging<br>Specifications | 500 mL, 5 L, 20 L |
|-----------------------------|-------------------|

|              |       |
|--------------|-------|
| Product Code | p1100 |
|--------------|-------|

## Intended Use

This product is a cleaning solution for the reaction system of the automated coagulation analyzer during the testing process, to facilitate the in vitro testing of the substance to be measured

## Principle

This product is a cleaning solution for the reaction system. Its main working principle is that the surfactant in the cleaning solution effectively lowers the surface tension of the proteins in the reaction system, so that the proteins in blood specimens are readily washed away...

This document is compiled from product registration information for reference only. Refer to the product instructions and the registration certificate for details.