

In Vitro Diagnostic Reagent

WF Cleaning Buffer Solution

An auxiliary reagent for immunoassay, used together with immunoassay reagents to provide/maintain the reaction environment for the Werfen series chemiluminescence immunoassay analyzers (systems)



Product Photo

Product Model

WF Cleaning Buffer Solution

Packaging Specifications

500 mL, 5 L, 10 L, 20 L

Product Code

p1118

Intended Use

An auxiliary reagent for immunoassay, used together with immunoassay reagents to provide/maintain the reaction environment for the Werfen series chemiluminescence immunoassay analyzers (systems)

Composition

Potassium chloride, phosphates, preservatives

Storage & Handling

Handle gently during loading and unloading, stack according to the arrow markings, and store in a dry, well-ventilated...

Other Information

Wash Buffer

