


In Vitro Diagnostic Reagent

HY Cleaning Buffer Solution

An auxiliary reagent for immunoassay, used together with immunoassay reagents to provide/maintain the reaction environment for the Changguang Huayi series chemiluminescence immunoassay analyzers



Product Photo

Product Model	HY Cleaning Buffer Solution
Packaging Specifications	500 mL, 500 mL × 2
Product Code	p1142

Intended Use

An auxiliary reagent for immunoassay, used together with immunoassay reagents to provide/maintain the reaction environment for the Changguang Huayi series chemiluminescence immunoassay analyzers

Composition

Composed of sodium chloride, phosphates and water, etc

Storage & Handling

Handle gently during loading and unloading, stack according to the arrow markings, and store in a dry, well-ventilated, clean indoor room at 2—...

This document is compiled from product registration information for reference only. Refer to the product instructions and the registration certificate for details.