

In Vitro Diagnostic Reagent

SL Concentrated Cleaning Buffer Solution

An auxiliary reagent for immunoassay, used together with immunoassay reagents to provide/maintain the reaction environment for the DiaSorin series chemiluminescence immunoassay analyzers (systems)



Product Photo

Product Model	SL Concentrated Cleaning Buffer Solution
Packaging Specifications	1 L
Product Code	p1152

Intended Use

An auxiliary reagent for immunoassay, used together with immunoassay reagents to provide/maintain the reaction environment for the DiaSorin series chemiluminescence immunoassay analyzers (systems)

Composition

Composed of buffering agent, preservative, surfactant and water

Storage & Handling

Handle gently during loading and unloading; stack according to the arrow markings; store in a dry, well-ventilated, clean indoor room at 2–35...

This document is compiled from product registration information for reference only. Refer to the product instructions and the registration certificate for details.