


In Vitro Diagnostic Reagent

WG Cleaning Buffer Solution

An auxiliary reagent for immunoassay, used together with immunoassay reagents to provide/maintain the reaction environment for the Weigao series chemiluminescence immunoassay analyzers



Product Photo

Product Model

WG Cleaning Buffer
SolutionPackaging
Specifications500 mL, 1 L, 5 L, 10 L, 20
L

Product Code

p1155

Intended Use

An auxiliary reagent for immunoassay, used together with immunoassay reagents to provide/maintain the reaction environment for the Weigao series chemiluminescence immunoassay analyzers

Composition

Composed of sodium chloride, phosphates, preservative

Storage & Handling

Handle gently during loading and unloading; stack according to the arrow markings; store in a dry, ventilated...

Other Information

Wash Buffer

