

In Vitro Diagnostic Reagent

Coagulation Analyzer Cleaning Solution–PLS Needle Wash 15 mL

This product is a cleaning solution for the reaction system during the testing process of the automated coagulation analyzer, so as to facilitate the in vitro testing of the substances to be measured



Product Photo

Product Model	Coagulation Analyzer PLS Needle Wash
Packaging Specifications	15 mL, 100 mL, 500 mL, 20 L
Product Code	p1225

Intended Use

This product is a cleaning solution for the reaction system during the testing process of the automated coagulation analyzer, so as to facilitate the in vitro testing of the substances to be measured

Principle

This product is a cleaning solution for the reaction system. Its main working principle is that the sodium hypochlorite in the cleaning solution can effectively decompose the proteins in the reaction system, making the proteins in the blood sample easier to clean...

This document is compiled from product registration information for reference only. Refer to the product instructions and the registration certificate for details.