

In Vitro Diagnostic Reagent

Coagulation Analyzer Cleaning Solution–BE Rinse

This product is a cleaning solution for the reaction system during the testing process of the automated coagulation analyzer, so as to facilitate the in vitro testing of the substances to be measured



Product Photo

Product Model	Coagulation Analyzer BE Rinse
Packaging Specifications	100 mL, 500 mL, 20 L
Product Code	p1226

Intended Use

This product is a cleaning solution for the reaction system during the testing process of the automated coagulation analyzer, so as to facilitate the in vitro testing of the substances to be measured

Principle

This product is a cleaning solution for the reaction system. Its main working principle is that the surfactants in the cleaning solution can effectively reduce the surface tension of the proteins in the reaction system...

This document is compiled from product registration information for reference only. Refer to the product instructions and the registration certificate for details.