

In Vitro Diagnostic Reagent

Coagulation Analyzer Cleaning Solution–BKN

This product is used to clean the reaction system during the testing process, so as to facilitate the in vitro testing of the substances to be measured



Product Photo

Product Model	Coagulation Analyzer Cleaning Solution–BKN
Packaging Specifications	500 mL, 3.3 L, 5 L, 20 L
Product Code	p1227

Intended Use

This product is used to clean the reaction system during the testing process, so as to facilitate the in vitro testing of the substances to be measured

Principle

This product is a system cleaning solution for the automated coagulation analyzer, which ensures efficient cleaning of the reagent needle and sample needle of the instrument and maintains the integrity of the hydraulic system, so as to ensure accurate aspiration volume and precise test results...

This document is compiled from product registration information for reference only. Refer to the product instructions and the registration certificate for details.