

In Vitro Diagnostic Reagent

Coagulation Analyzer Cleaning Solution–JS

This product is a cleaning solution for the reaction system during the testing process of the automated coagulation analyzer, so as to facilitate the in vitro testing of the substances to be measured



Product Photo

Product Model	Coagulation Analyzer Cleaning Solution–JS
Packaging Specifications	250 mL, 500 mL, 1 L, 20 L
Product Code	p1228

Intended Use

This product is a cleaning solution for the reaction system during the testing process of the automated coagulation analyzer, so as to facilitate the in vitro testing of the substances to be measured

Principle

This product is a cleaning solution for the reaction system. Its main working principle is that the potassium hydroxide in the cleaning solution can dissolve the proteins in the reaction system, making the proteins in the blood sample easier to clean...

This document is compiled from product registration information for reference only. Refer to the product instructions and the registration certificate for details.