

In Vitro Diagnostic Reagent

Sample Diluent-AU

Used to dilute and liquefy the sample to be tested, so that the substances to be measured can be assayed with the electrolyte module. The reagent itself does not directly participate in the measurement.



Product Photo

Product Model	Sample Diluent-AU
Packaging Specifications	2.0 L, 2 L x 2
Product Code	68

Intended Use

Used to dilute and liquefy the sample to be tested, so that the substances to be measured can be assayed with the electrolyte module. The reagent itself does not directly participate in the measurement.

Compatible Devices

Beckman Coulter AU5800

This document is compiled from product registration information for reference only. Refer to the product instructions and the registration certificate for details.