

In Vitro Diagnostic Reagent

Sample Diluent-RL

Used to dilute and liquefy the sample to be tested, so that the substances to be measured can be assayed with the electrolyte module. The reagent itself does not directly participate in the measurement.



Product Photo

Product Model	Sample Diluent-RL
Packaging Specifications	2.0 L, 2 L × 2, 300 mL
Product Code	72

Intended Use

Used to dilute and liquefy the sample to be tested, so that the substances to be measured can be assayed with the electrolyte module. The reagent itself does not directly participate in the measurement.

Compatible Devices

- Hitachi Hitachi 3100
- Hitachi Hitachi 7020
- Hitachi Hitachi 7060
- Hitachi Hitachi 7080
- Hitachi Hitachi 7100
- Hitachi Hitachi 7150
- Hitachi Hitachi 7170
- Hitachi Hitachi 7180
- Hitachi Hitachi 7600
- Hitachi LABOSPECT 008

This document is compiled from product registration information for reference only. Refer to the product instructions and the registration certificate for details.