

In Vitro Diagnostic Reagent

Electrode Cleaning Solution–RL 100 mL

Used to dilute and liquefy the sample to be tested, so that the substances to be measured can be assayed with the electrolyte module. The reagent itself does not directly participate in the measurement



Product Photo

Product Model	Sample Diluent–RH
Packaging Specifications	2.0 L, 2 Lx2, 300 mL
Product Code	p1064

Intended Use

Used to dilute and liquefy the sample to be tested, so that the substances to be measured can be assayed with the electrolyte module. The reagent itself does not directly participate in the measurement

Principle

The electrode method for determining electrolytes is based on the Nern...

This document is compiled from product registration information for reference only. Refer to the product instructions and the registration certificate for details.