

In Vitro Diagnostic Reagent

Hemolysing Agent-BC II

Used prior to blood cell analysis to lyse red blood cells, release hemoglobin and preserve the morphology of the cells to be analyzed, thus facilitating cell differentiation and counting or quantitative determination of hemoglobin.



Product Photo

Product Model	Hemolysing Agent-BC II
Packaging Specifications	500 mL, 1 L, 5 L, 20 L
Product Code	146

Intended Use

Used prior to blood cell analysis to lyse red blood cells, release hemoglobin and preserve the morphology of the cells to be analyzed, thus facilitating cell differentiation and counting or quantitative determination of hemoglobin.

Compatible Devices

Mindray BC-1800

Mindray BC-2300

Mindray BC-2600

Mindray BC-2800

Mindray BC-2900

Mindray BC-3000 Plus

Mindray BC-3200

Hualiang HA series

Meisikang MC series

Pukang PE-6000

Pukang PE-6100

Pukang PE-6300

Pukang PE-6800

Pukang PE-6900

