

In Vitro Diagnostic Reagent

Hemolysing Agent–EO II 53

Used prior to blood cell analysis to lyse red blood cells, release hemoglobin and preserve the morphology of the cells to be analyzed, thus facilitating cell differentiation and counting or quantitative determination of hemoglobin.



Product Photo

Product Model	Hemolysing Agent–EO II 53
Packaging Specifications	200 mL, 500 mL, 1 L, 5 L, 20 L
Product Code	152

Intended Use

Used prior to blood cell analysis to lyse red blood cells, release hemoglobin and preserve the morphology of the cells to be analyzed, thus facilitating cell differentiation and counting or quantitative determination of hemoglobin.

Compatible Devices

Mindray BC–5100

Mindray BC–5180

Mindray BC–5300

Mindray BC–5310

Mindray BC–5380

This document is compiled from product registration information for reference only. Refer to the product instructions and the registration certificate for details.

