

In Vitro Diagnostic Reagent

Hemolysing Agent–EO II 53S

Used prior to blood cell analysis to lyse red blood cells, release hemoglobin and preserve the morphology of the cells to be analyzed, thus facilitating cell differentiation and counting or quantitative determination of hemoglobin.



Product Photo

|                          |                                |
|--------------------------|--------------------------------|
| Product Model            | Hemolysing Agent–EO II 53S     |
| Packaging Specifications | 200 mL, 500 mL, 1 L, 5 L, 20 L |
| Product Code             | 156                            |

Intended Use

Used prior to blood cell analysis to lyse red blood cells, release hemoglobin and preserve the morphology of the cells to be analyzed, thus facilitating cell differentiation and counting or quantitative determination of hemoglobin.

Compatible Devices

- Mindray BC–5100
- Mindray BC–5180
- Mindray BC–5380
- Mindray BC–5390

This document is compiled from product registration information for reference only. Refer to the product instructions and the registration certificate for details.