

In Vitro Diagnostic Reagent

Hemolysing Agent-5DIFF

Used prior to blood cell analysis to lyse red blood cells, release hemoglobin and preserve the morphology of the cells to be analyzed; when required, the white blood cells are stained before cell differentiation and counting, or cell differentiation and counting and/or quantitative determination of hemoglobin are carried out directly, etc. ...



Product Photo

Product Model	Hemolysing Agent-5DIFF
Packaging Specifications	Hemolysing Agent Pack I 1.9 L + Hemolysing Agent Pack II 850 mL
Product Code	p1109

Intended Use

Used prior to blood cell analysis to lyse red blood cells, release hemoglobin and preserve the morphology of the cells to be analyzed; when required, the white blood cells are stained before cell differentiation and counting, or cell differentiation and counting and/or quantitative determination of hemoglobin are carried out directly, etc. ...

This document is compiled from product registration information for reference only. Refer to the product instructions and the registration certificate for details.