

In Vitro Diagnostic Reagent

Hemolysing Agent-LD2 500 mL

Used prior to blood cell analysis to lyse red blood cells, release hemoglobin and preserve the morphology of the cells to be analyzed; when required, the white blood cells are stained before cell differentiation and counting, or cell differentiation and counting and/or quantitative determination of hemoglobin are carried out directly



Product Photo

Product Model	Hemolysing Agent-LD2
Packaging Specifications	100 mL, 200 mL, 500 mL, 1 L, 4 L, 10 L, 20 L
Product Code	p1127

Intended Use

Used prior to blood cell analysis to lyse red blood cells, release hemoglobin and preserve the morphology of the cells to be analyzed; when required, the white blood cells are stained before cell differentiation and counting, or cell differentiation and counting and/or quantitative determination of hemoglobin are carried out directly

This document is compiled from product registration information for reference only. Refer to the product instructions and the registration certificate for details.