

In Vitro Diagnostic Reagent

Hemolysing Agent-LR3

Used prior to blood cell analysis to lyse red blood cells, release hemoglobin and preserve the morphology of the cells to be analyzed; when required, the white blood cells are stained before cell differentiation and counting, or cell differentiation and counting and/or quantitative determination of hemoglobin are carried out directly



Product Photo

Product Model	Hemolysing Agent-LR3
Packaging Specifications	200 mL, 500 mL, 1 L, 4 L, 5 L, 10 L, 20 L
Product Code	p1129

Intended Use

Used prior to blood cell analysis to lyse red blood cells, release hemoglobin and preserve the morphology of the cells to be analyzed; when required, the white blood cells are stained before cell differentiation and counting, or cell differentiation and counting and/or quantitative determination of hemoglobin are carried out directly

Principle

It lyses red blood cells while keeping the morphology of white blood cells...

This document is compiled from product registration information for reference only. Refer to the product instructions and the registration certificate for details.

