

In Vitro Diagnostic Reagent

Hemolysing Agent-UT II

Used prior to blood cell analysis to lyse red blood cells, release hemoglobin and preserve the morphology of the cells to be analyzed, thus facilitating cell differentiation and counting or quantitative determination of hemoglobin



Product Photo

Product Model	Hemolysing Agent-BC II
Packaging Specifications	500 mL, 1 L, 5 L, 20 L
Product Code	p1130

Intended Use

Used prior to blood cell analysis to lyse red blood cells, release hemoglobin and preserve the morphology of the cells to be analyzed, thus facilitating cell differentiation and counting or quantitative determination of hemoglobin

Principle

It lyses red blood cells while keeping the morphology of white blood cells unchanged, enabling the analyzer to count and classify white blood cells by the impedance method

Composition

Quaternary ammonium salt, buffering agent...

