

In Vitro Diagnostic Reagent

Sheath Fluid–UT5S 20 L

Used together with a cell analyzer; it dilutes the blood sample to form a sheath stream, facilitating cell counting and classification by the analyzer



Product Photo

Intended Use

Used together with a cell analyzer; it dilutes the blood sample to form a sheath stream, facilitating cell counting and classification by the analyzer

Principle

It dilutes the white blood cells, enabling the analyzer to count and classify white blood cells by the laser light scatter method

Composition

Hydrochloric acid, tris(hydroxymethyl)aminomethane...

Product Model	Sheath Fluid–UT5S
Packaging Specifications	500 mL, 1 L, 5 L, 10 L, 20 L
Product Code	p1132

This document is compiled from product registration information for reference only. Refer to the product instructions and the registration certificate for details.