

In Vitro Diagnostic Reagent

Diluent-RD

Used to dilute the whole blood sample to be tested and, together with the blood analyzer or the reagent, to test the analyte. Clinically, it is used prior to blood cell analysis to prepare a cell suspension by diluting the sample



Product Photo

Product Model

Diluent-RD

Packaging
Specifications

1 L, 2 L, 5 L, 10 L, 20 L

Product Code

p1138

Intended Use

Used to dilute the whole blood sample to be tested and, together with the blood analyzer or the reagent, to test the analyte. Clinically, it is used prior to blood cell analysis to prepare a cell suspension by diluting the sample

Principle

It dilutes the blood cells, enabling the analyzer to count and classify the blood cells by the impedance method

Composition

Sodium chloride, boric acid...

