

In Vitro Diagnostic Reagent

Stain Solution–MDF II

Used to stain blood cells so that their morphology and structure can be observed, enabling the blood analyzer to perform blood cell differentiation and counting



Product Photo

Product Model	Stain Solution–MDF II
Packaging Specifications	12 mL, 22 mL, 22 mLx2, 42 mL, 42 mLx2, 1 L, 2 L
Product Code	p1188

Intended Use

Used to stain blood cells so that their morphology and structure can be observed, enabling the blood analyzer to perform blood cell differentiation and counting

Principle

White blood cells are stained with the aid of the stain solution for blood cell analysis and, in combination with the flow cytometric technique based on laser light scatter and fluorescence...

This document is compiled from product registration information for reference only. Refer to the product instructions and the registration certificate for details.