

In Vitro Diagnostic Reagent

Sheath Fluid-CD

For use with cell analyzers; it is used to dilute the blood sample to form a sheath flow, facilitating cell counting and classification by the analyzer



Product Photo

Product Model	Sheath Fluid-CD
Packaging Specifications	10 L, 20 L
Product Code	p1189

Intended Use

For use with cell analyzers; it is used to dilute the blood sample to form a sheath flow, facilitating cell counting and classification by the analyzer

Principle

It dilutes the white blood cells, enabling the analyzer to count and classify the white blood cells by the laser light scatter method...

Other Information

Sheath fluid

This document is compiled from product registration information for reference only. Refer to the product instructions and the registration certificate for details.