

## In Vitro Diagnostic Reagent

## Hemolysing Agent-M

Used prior to blood cell analysis to lyse red blood cells, release hemoglobin and preserve the morphology of the cells to be analyzed, thus facilitating cell differentiation and counting or quantitative determination of hemoglobin



Product Photo

|                          |                    |
|--------------------------|--------------------|
| Product Model            | Hemolysing Agent-M |
| Packaging Specifications | 5 L, 20 L          |
| Product Code             | p1201              |

## Intended Use

Used prior to blood cell analysis to lyse red blood cells, release hemoglobin and preserve the morphology of the cells to be analyzed, thus facilitating cell differentiation and counting or quantitative determination of hemoglobin

## Principle

It lyses red blood cells while keeping the morphology of white blood cells unchanged, enabling the analyzer to count and classify white blood cells by the impedance method

## Composition

Quaternary ammonium salt, buffering agent...

