

In Vitro Diagnostic Reagent

Diluent-KT

Used to dilute the sample and prepare a cell suspension prior to blood cell analysis



Product Photo

Product Model	Diluent-KT
Packaging Specifications	5 L, 10 L, 20 L
Product Code	p1205

Intended Use

Used to dilute the sample and prepare a cell suspension prior to blood cell analysis

Principle

It dilutes the blood cells, enabling the analyzer to count and classify the blood cells by the impedance method

Composition

Sodium sulfate, phosphates, preservatives

Storage & Handling

Handle with care during loading and unloading; stack according to the arrow markings; store in a dry, well-ventilated, clean indoor room at 2...

