

In Vitro Diagnostic Reagent

Diluent-KL2

Used to dilute the whole blood sample to be tested and, together with the blood analyzer or the reagent, to test the analyte. Clinically, it is used prior to blood cell analysis to prepare a cell suspension by diluting the sample



Product Photo

Product Model

Diluent-KL2

Packaging Specifications

1 L, 2 L, 5 L, 10 L, 20 L

Product Code

p1209

Intended Use

Used to dilute the whole blood sample to be tested and, together with the blood analyzer or the reagent, to test the analyte. Clinically, it is used prior to blood cell analysis to prepare a cell suspension by diluting the sample

Principle

It dilutes the blood cells, enabling the analyzer to count and classify the blood cells by the impedance method

Composition

Sodium hydroxide, sodium sulfate, disodium ethylenediaminetetraacetate...

