

In Vitro Diagnostic Reagent

Hemolysing Agent-90EOS

Used prior to blood cell analysis to lyse red blood cells, release hemoglobin and preserve the morphology of the cells to be analyzed, thus facilitating cell differentiation and counting or quantitative determination of hemoglobin



Product Photo

Product Model	Hemolysing Agent-90EOS
Packaging Specifications	500 mL, 1 L × 2, 5 L, 10 L, 20 L
Product Code	p1215

Intended Use

Used prior to blood cell analysis to lyse red blood cells, release hemoglobin and preserve the morphology of the cells to be analyzed, thus facilitating cell differentiation and counting or quantitative determination of hemoglobin

Principle

During blood cell testing, this reagent lyses red blood cells and stains the white blood cells; after a period of time, the staining reaction will be...

This document is compiled from product registration information for reference only. Refer to the product instructions and the registration certificate for details.