

In Vitro Diagnostic Reagent

Hemolysing Agent-Z5LB 500 mL

Used prior to blood cell analysis to lyse red blood cells, release hemoglobin and preserve the morphology of the cells to be analyzed; when required, the white blood cells are stained before cell differentiation and counting, or cell differentiation and counting and/or quantitative determination of hemoglobin are carried out directly



Product Photo

Product Model	Hemolysing Agent-Z5LB
Packaging Specifications	50 mL, 100 mL, 100 mL × 4, 200 mL, 500 mL, 1 L, 5 L, 10 L, 20 L
Product Code	p1217

Intended Use

Used prior to blood cell analysis to lyse red blood cells, release hemoglobin and preserve the morphology of the cells to be analyzed; when required, the white blood cells are stained before cell differentiation and counting, or cell differentiation and counting and/or quantitative determination of hemoglobin are carried out directly

Principle

It lyses red blood cells...

This document is compiled from product registration information for reference only. Refer to the product instructions and the registration certificate for details.

