

In Vitro Diagnostic Reagent

Hemolysing Agent-Z3

Used prior to blood cell analysis to lyse red blood cells, release hemoglobin and preserve the morphology of the cells to be analyzed; when required, the white blood cells are stained before cell differentiation and counting, or cell differentiation and counting and/or quantitative determination of hemoglobin are carried out directly



Product Photo

Product Model	Hemolysing Agent-Z3
Packaging Specifications	75 mL, 150 mL, 250 mL, 500 mL, 1 L, 5 L, 10 L, 20 L
Product Code	p1222

Intended Use

Used prior to blood cell analysis to lyse red blood cells, release hemoglobin and preserve the morphology of the cells to be analyzed; when required, the white blood cells are stained before cell differentiation and counting, or cell differentiation and counting and/or quantitative determination of hemoglobin are carried out directly

Principle

It lyses red blood cells while keeping the white blood cells...

This document is compiled from product registration information for reference only. Refer to the product instructions and the registration certificate for details.

