

In Vitro Diagnostic Reagent

Diluent-DM

Used to dilute the sample and prepare a cell suspension prior to blood cell analysis



Product Photo

Product Model	Diluent-DM
Packaging Specifications	5 L, 20 L
Product Code	p1266

Intended Use

Used to dilute the sample and prepare a cell suspension prior to blood cell analysis

Principle

It dilutes the blood cells, enabling the analyzer to count and classify the blood cells by the impedance method

Composition

Sodium chloride, sodium sulfate, disodium ethylenediaminetetraacetate, preservatives

Storage & Handling

Handle with care during loading and unloading; stack according to the arrow markings; store in a dry, well-ventilated...

