

In Vitro Diagnostic Reagent

Diluent-FDL

Used to dilute the sample and prepare a cell suspension prior to blood cell analysis



Product Photo

|                          |                                  |
|--------------------------|----------------------------------|
| Product Model            | Diluent-FDL                      |
| Packaging Specifications | 1 L, 1.5 L, 4 L, 5 L, 10 L, 20 L |
| Product Code             | p1285                            |

Intended Use

Used to dilute the sample and prepare a cell suspension prior to blood cell analysis

Principle

It dilutes the blood cells, enabling the analyzer to count and classify the blood cells by the impedance method

Composition

Buffering agent, surfactant

Storage & Handling

Handle with care during loading and unloading; stack according to the arrow markings...

This document is compiled from product registration information for reference only. Refer to the product instructions and the registration certificate for details.