

## In Vitro Diagnostic Reagent

## Hemolysing Agent-XKFB

Used prior to blood cell analysis to lyse red blood cells, release hemoglobin and preserve the morphology of the cells to be analyzed, thus facilitating cell differentiation and counting or quantitative determination of hemoglobin



Product Photo

|                          |                       |
|--------------------------|-----------------------|
| Product Model            | Hemolysing Agent-XKFB |
| Packaging Specifications | 500 mL, 5 L, 20 L     |
| Product Code             | p1288                 |

## Intended Use

Used prior to blood cell analysis to lyse red blood cells, release hemoglobin and preserve the morphology of the cells to be analyzed, thus facilitating cell differentiation and counting or quantitative determination of hemoglobin

## Principle

It lyses red blood cells while keeping the morphology of white blood cells unchanged, enabling the analyzer to count and classify white blood cells by the laser light-scattering method

This document is compiled from product registration information for reference only. Refer to the product instructions and the registration certificate for details.