


In Vitro Diagnostic Reagent

Hemolysing Agent–MDF

Used prior to blood cell analysis to lyse red blood cells, release hemoglobin and preserve the morphology of the cells to be analyzed, thus facilitating cell differentiation and counting or quantitative determination of hemoglobin



Product Photo

Product Model

Hemolysing Agent–MDF

Packaging Specifications

500 mL, 4 L, 5 L, 20 L

Product Code

p1289

Intended Use

Used prior to blood cell analysis to lyse red blood cells, release hemoglobin and preserve the morphology of the cells to be analyzed, thus facilitating cell differentiation and counting or quantitative determination of hemoglobin

Principle

It lyses red blood cells while keeping the morphology of white blood cells unchanged, enabling the analyzer to count and classify white blood cells by the laser light–scattering method

This document is compiled from product registration information for reference only. Refer to the product instructions and the registration certificate for details.