

In Vitro Diagnostic Reagent

Stain Solution–MDF

Used to stain blood cells so that their morphology and structure can be observed, enabling the blood analyzer to perform blood cell differentiation and counting



Product Photo

Product Model	Stain Solution–MDF
Packaging Specifications	22 mL, 42 mL, 82 mL
Product Code	p1293

Intended Use

Used to stain blood cells so that their morphology and structure can be observed, enabling the blood analyzer to perform blood cell differentiation and counting

Principle

White blood cells are stained with the aid of the stain solution for blood cell analysis, and in combination with the flow cytometric technique based on laser light scatter and fluorescence staining, white blood cell detection is performed...

This document is compiled from product registration information for reference only. Refer to the product instructions and the registration certificate for details.