

In Vitro Diagnostic Reagent

Stain Solution–FN68P

Used to stain blood cells so that their morphology and structure can be observed, enabling the blood analyzer to perform blood cell differentiation and counting



Product Photo

|                          |                                              |
|--------------------------|----------------------------------------------|
| Product Model            | Stain Solution–FN68P                         |
| Packaging Specifications | 12 mL, 12 mL × 4, 48 mL, 48 mL × 4, 1 L, 2 L |
| Product Code             | p1303                                        |

Intended Use

Used to stain blood cells so that their morphology and structure can be observed, enabling the blood analyzer to perform blood cell differentiation and counting

Principle

During blood cell testing, Stain Solution–FN68P, with the aid of Hemolysing Agent–LN68P, acts on the nucleated red blood cells and white blood cells...

This document is compiled from product registration information for reference only. Refer to the product instructions and the registration certificate for details.