

In Vitro Diagnostic Reagent

Hemolysing Agent-MB-5EO II

Used prior to blood cell analysis to lyse red blood cells, release hemoglobin and preserve the morphology of the cells to be analyzed; when required, the white blood cells are stained before cell differentiation and counting, or cell differentiation and counting and/or quantitative determination of hemoglobin are carried out directly



Product Photo

Product Model	Hemolysing Agent-MB-5EO II
Packaging Specifications	500 mL, 1 L, 5 L, 10 L, 20 L
Product Code	p1315

Intended Use

Used prior to blood cell analysis to lyse red blood cells, release hemoglobin and preserve the morphology of the cells to be analyzed; when required, the white blood cells are stained before cell differentiation and counting, or cell differentiation and counting and/or quantitative determination of hemoglobin are carried out directly

Principle

It lyses red blood cells while keeping the morphology of white blood cells unchanged...

This document is compiled from product registration information for reference only. Refer to the product instructions and the registration certificate for details.

