

In Vitro Diagnostic Reagent

Sheath Fluid–NUT

For use with analyzers that employ the principle of flow cytometric analysis. It is used to dilute urine samples to form a sheath flow, facilitating cell counting and classification by the analyzer.



Product Photo

Product Model	Sheath Fluid–NUT
Packaging Specifications	500 mL, 1 L, 5 L, 10 L, 15 L, 20 L
Product Code	712

Intended Use

For use with analyzers that employ the principle of flow cytometric analysis. It is used to dilute urine samples to form a sheath flow, facilitating cell counting and classification by the analyzer.

Compatible Devices

Urit USS–80

This document is compiled from product registration information for reference only. Refer to the product instructions and the registration certificate for details.