

In Vitro Diagnostic Reagent

Sheath Fluid-NUF

For use with analyzers that employ the principle of flow cytometric analysis. It is used to dilute urine samples to form a sheath flow, facilitating cell counting and classification by the analyzer.



Product Photo

Product Model	Sheath Fluid-NUF
Packaging Specifications	500 mL, 1 L, 10 L, 20 L
Product Code	715

Intended Use

For use with analyzers that employ the principle of flow cytometric analysis. It is used to dilute urine samples to form a sheath flow, facilitating cell counting and classification by the analyzer.

Compatible Devices

- Sysmex UF-3000
- Sysmex UF-4000
- Sysmex UF-5000

This document is compiled from product registration information for reference only. Refer to the product instructions and the registration certificate for details.